



PathoPlus Express - Getting Started Guide

Installation and Setup


1. Download Setup of PathoPlus Express from www.charlycomputers.com. Extract the downloaded Zip file by archival programs like WinZip or WinRAR.
2. PathoPlusExpress_Setup.exe will be extracted from the Zip file. Double-click and run the setup file.
3. If a User Access Control (UAC) confirmation dialogue appears, click “Yes” and allow setup.exe to execute. The following screen will appear. Click on “OK”.




4. Click on the button shown in figure to continue installation.




5. Click on Continue button.




6. Program will start installing and progress bar will be shown.




7. If you are installing PathoPlus Express on Windows 7 and above, you may encounter Version Conflict dialog asking to keep the newer file as shown in figure below. You may encounter 4 to 5 such dialogs depending upon your systems configuration. Click on Yes button for all such dialogs and your installation will continue.




8. If you are installing PathoPlus Express on Windows 7 and above, you may encounter a message as shown in below figure. Click on OK.




Then in the below screen click on Ignore.




9. PathoPlus Express Setup was completed successfully message will be displayed after completion of setup.



10. Keep Internet connected and run the program from Start > All Programs > PathoPlus > PathoPlus Express
11. Your product copy will be activated instantly for a trial period of 30 days. Click on Exit.



12. PathoPlus Express will launch and you may disconnect from internet if you wish.
Click on Ok and your program will start.



Purchasing License


1. You can purchase license for PathoPlus Express either before or after the evaluation/demo period expires.
2. Click on "Buy Now" button on PathoPlus Express page to make payment for your license.
3. Enter your details (Name, Email, Contact Number etc.) and your ClientID on the Buy page.
4. Make payment through Credit Card/Debit Card/Net Banking.
5. After successful payment, your ClientID will be instantly activated.
6. Fill registration form in PathoPlus Express and click on "Register" button to complete the license activation process.

Note: - Internet connection will be required only at first run and license registration. License validity will be lifetime with single user on single computer.

How to find My Client ID?

If you have not downloaded PathoPlus Express Demo, then first download and install PathoPlus Express Demo on your computer.

- 1) If your Demo has not expired then you can find your Client ID on Login screen as shown below.



- 2) If your Demo has expired then you can find your Client ID on the message screen as shown below.

